

How to use “FootScan” Foot Scanner

The Scanning of a client’s foot is done in “Three parts”.

1st part is to take image of the foot, 2nd part is to fill client’s data & 3rd part is to upload data of client to our website.

Part I

Start by connecting the USB of the scanner with your Desktop /Laptop computer.

Double Click FootScan button on Desktop. The below Screen will open in which a small window will appear “click connect button”, click **OK**

The screenshot shows the FootScan software interface. The main window has a menu bar with 'File', 'Connect', 'Pre-Scan', 'Scan', 'Edit Client', 'Client Data', and 'Help'. A dialog box titled 'Click Connect Button.' with an 'OK' button is centered on the screen. The background interface includes a large empty area for the foot scan, a 'Foot Pressure %' section with input fields for Forefoot, Heel, Left Foot, and Right Foot, and a section for foot measurements (Arch, Arch Type, Length, Width, Foot Type, Shoe Size, Ball Joint) for both Left and Right feet. At the bottom, there are fields for Date, Name, Last Name, Sex, I.D., SR No, Mobile, Email, Country, State, City, Pin, Address, Tel Code, DOB, Height, Weight, Diabetic, and Recommended Model, along with a 'Save' button. A 'Save All' button is also present at the bottom right. A note at the bottom right reads: 'Note: Print this foot scan report for consultation with a certified doctor/podiatrist'.


The top left “Connect” button will become enabled by turning “Pink”. Click this button.

The screenshot shows the FootScan software interface with the 'Connect' button in the menu bar highlighted with a red circle. The rest of the interface is the same as in the previous screenshot, including the main scan area, measurement sections, and client data entry fields.

A small window opens




After a while a small window appears “click Pre- scan button”. Click **OK**




The “**Pre-scan**” button will become enabled by turning “Pink”. Click this button & a message appears



After a moment a small window opens “Please Stand bare foot on the machine”. ASK YOUR CLIENT TO STAND BARE FOOT ON THE SCANNER in an Upright posture & Click **OK**



A small window appears “Click scan button”, Click




On the top the “Scan” button will become enabled by turning “Pink”. Click this button & in a few seconds a “Foot Image of the client” appears on the screen along with foot details & a small window opens “scanning complete”.

If you or your client are satisfied with the result, click & proceed to fill client DATA in the boxes given below.

“BUT” if you or your client feels that the image appearing on the screen is not correct, then ask your client to remain standing on the scanner & click button.

The top “Re-scan” will become enabled by becoming “Pink”. Click this button to again scan foot of your client once again & follow earlier process.



The screenshot shows the software interface with a menu bar (File, Disconnect, Pre-Scan, Re-Scan, Edit Client, Client Data, Help) and a main window. The main window displays two foot pressure heatmaps (left and right) and a camera view of the scanner. A 'Confirm' dialog box is open, asking for confirmation to proceed. The right panel shows 'Foot Pressure %' data for Left and Right feet, including Forefoot, Heel, Arch, and Ball Joint measurements. The bottom panel contains client data entry fields for Date, Name, Last Name, Sex, I.D., SR No, Mobile, Email, Country, State, City, Pin, Address, Tel Code, DOB, Height, Weight, Diabetic status, and Recommended Model. A 'Save All' button is located at the bottom right.

Foot Pressure %	
Forefoot	55.85 %
Heel	44.15 %
Left Foot	58.52 %
Right Foot	41.48 %
Arch	50.59 %
Ball Joint	0.04 %
Forefoot	83.05 %
Heel	16.95 %

Arch Type: Left Foot: Low, Right Foot: High

Foot Type: Left Foot: Select, Right Foot: Select

Shoe Size: Left Foot: [], Right Foot: []

Ball Joint: Left Foot: [] mm, Right Foot: [] mm (Lateral/Medial)


Client Data: Date: 06-Jun-2019, Name: [], Last Name: [], Sex: [], I.D.: 2, SR No: [], Mobile: [], Email: [], Country: [], State: [], City: [], Pin: [], Address: [], Tel Code: [], DOB: [], Height: [] cm, Weight: [] kg, Diabetic: [], Recommended Model: []

Note: Print this foot scan report for consultation with a certified doctor/podiatrist


Part II

Click **YES** in the small window “scanning complete” & proceed below to fill client’s personal details

After filling all personal details in the boxes in the bottom side of window, click **SAVE** button. A small window opens “**Basic detail of client saved successfully**”. This window gives **YES** **NO** option. Click **YES** to proceed filling all other details of foot in the right side boxes & click **SAVE ALL** button to complete a person’s foot scanning.



A small window appears stating Data Stored successfully. Click **OK** to complete this foot scanning.



BUT

for any reason if you want to change client’s Personal details, then click **NO** & go back to **Edit client** button on top & change clients personal & other foot details & follow all other procedure as mentioned earlier.

Part III

Click “client data” button on top right side of the screen to open “Upload scanning” page. In this page you will find all foot scanning reports of people. Here you can **Edit/print/delete** details of all clients before uploading it to our website.

On the left side of this page you’ll find small check boxes. Tick the check box **select all** & click **upload** to upload all scanning reports to us in one click. You can also select them one by one & upload them.

Now a small window will open “Data upload successful”. Click **OK**. All tick marked reports will shift to lower side of the center line indicating successful uploading.

This completes the entire scanning process & you can close the application now.

The screenshot shows the 'Upload Scannings' application window. At the top, there are search filters for 'Scan Date From' (06-Jun-2019) and 'To' (06-Jun-2019), with 'Search' and 'Clear' buttons. There are also checkboxes for 'Search in Pending' and 'Search in Uploaded', and buttons for 'Search Customer', 'Upload', and 'Export to Excel'. Below the filters is a table with the following data:

Select	Scan Date	SR No	Customer Name	City	Email	Mobile	Scan No	Action	Action	Action
<input checked="" type="checkbox"/>	06-Jun-19	2	Sarjay Singh Rana	Chandausi	ranasanjay111@gmail.com	9718381691	00002061000001	Edit	Print	Delete

A dialog box titled 'Data Uploaded Successful' is centered on the screen, with an 'OK' button. The background of the application window is dimmed. On the right side of the window, there are some percentage indicators and labels like 'Lateral' and 'Medial'. At the bottom of the window, there is a note: 'Note: Print this foot scan report for consultation with a certified doctor/podiatrist'.